HYDROS

WATER TREATMENT SYSTEM

Handbuch

Manual

Mode d’emploi

Handboek

Manuale

Manual

Podręcznik

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1 Definition of Terms

- ① Pressure Vessel
- ② Filter Cartridge
- ③ Pressure Vessel Lid
- ④ Connector Head
- ⑤ Inlet Hose
- ⑥ Inlet Valve on the Inlet Hose
- ⑦ Outlet Hose to the End Device
- ⑧ Locking
- ⑨ Kick Loop
- ⑩ Ejector Base
- ⑪ Mantle Handle
- ⑫ Transport Protective Cap
- ⑬ Lid Handle
- ⑭ Reducer 1"-3/4"
- ⑮ Flush valve with flush hose

2 General Information

2.1 Function and Application

The HYDROS Treatment System optimises tap water by means of full demineralisation.

In the HYDROS Treatment System calcium and magnesium ions are selectively removed from the drinking water via an ion exchanger using the flow method. This prevents residues on the dental instruments during drying due to carbonate hardness and scale deposits on important machine parts.

Furthermore, the HYDROS Treatment System removes salts in the water by means of full demineralisation of the tap water, preventing them from being deposited on instruments.

2.2 Guarantee Provisions

The HYDROS Treatment System is subject to the statutory guarantee of 2 years. A guarantee claim may be asserted only if all instructions in this Handbook are followed and observed.

2.3 Storage/Transport

Adhere to the ambient conditions in the Technical Data (Chapter10) for storage and transport.

The handbook should be seen as part of the product and kept for the whole service life of the filter system and passed on to subsequent owners.

2.4 Recycling/Disposal

By disposing of this product and its packaging in the correct manner you are helping to prevent potential negative impacts on people and the environment that could be caused by incorrect disposal. Comply with local regulations and bring the units to be disposed of to a proper recycling facility.

3 Operating and Safety Information

3.1 Qualified Personnel

Installation and maintenance of the filter system may be carried out only by trained or authorised personnel.

3.2 Correct Use

The product can only be operated properly and safely if it is installed, used and serviced in the manner described in this Manual. Use only filter cartridges suitable for this appliance.

Note: the HYDROS Treatment System and system-specific filter cartridges may only be used with the C61wd G4 and M2 G4 units.

3.3 Liability Exclusion

Installation must be carried out in accordance with the instructions in this Manual. SciCan shall not be held liable for any damage, including subsequent damage, arising from the incorrect installation or use of the product.

3.4 Specific Safety Information

- Only drinking water within the intake water temperature range listed in Chapter 10 may be used as intake water for the HYDROS water filter system. No microbiologically impaired water or water of unknown microbiological quality may be used without appropriate sterilisation.
- SciCan recommends that the filter system shall not be decommissioned for a long period. After stagnation times of over 4 weeks, the tank should be flushed or a pre-rinse executed. Please also note the maximum usage period of the filter cartridge is 12 months (Chapter 6).
- The filter system is not resistant to heavily concentrated cleaning agents (e.g. bleach, chlorinated solvents, heavy oxidants) and must not come into contact with them.
- The filter system must not be opened or dismantled during operation. The filter cartridge must not be opened.
- The pressure vessel and the pressure vessel lid of the filter systems have a service life of up to ten years (from the date of installation), provided that they are installed and used correctly and the operating conditions outlined in the Technical Data chapter are adhered to. They must always be replaced after a maximum of ten years. The hoses must be replaced in rotation after a maximum of five years.

- **Production date:**

Production code sticker filter cartridge and packaging - Example: B715002010

7	Production year, here: 2007
15	Production week, here: calendar week 15
002	Batch No. filter medium, here the second batch filled in terms of quantity
010	Serial number of the filter cartridge, here the tenth cartridge from the second batch

Production code sticker connector head - Example: 7102 XX

7	Production year, here: 2007
10	Production week, here: calendar week 10
2	Production day from Monday (1) to Friday (5), here: Tuesday
xx	Internal code

Production date pressure vessel and pressure vessel lid - Example: 0307

03	Production month, here: March
07	Production year, here: 2007

3.5 Technical Safety Assembly Instructions

- The terminal device operated with the filter must be free of limescale prior to installation. **⚠ Caution:** No water softening system may be operated upstream from the filter system.
- Protect the filter system from sunlight and mechanical damage. Do not assemble near sources of heat and open flames.
- A stop valve must be installed before the filter system intake hose.
- If the water pressure is higher than 6 bar, a pressure reducer must be installed before the filter system.
- A non-return valve tested by the DVGW has been factory-installed at the water intake of the filter head.
- No copper pipes and no galvanised or nickel-plated pipes/connectors may be installed between the water filter and the consumer. When choosing the material for parts that come into contact with water after the HYDROS filter system it must be remembered that, due to the process, decarbonised and/or fully demineralised water contains free carbon dioxide. The use of the HYDROS hose sets is recommended.
- All parts must be installed in accordance with the country-specific guidelines on the installation of water facilities.

4 Installation

⚠ Caution: Prior to installation read the technical data (Chapter 10) and the operating and safety information (Chapter 3). After the product has been stored at temperatures below 0°C, it must be stored (with the original packaging open) at the ambient temperature of the place of installation for at least 24 hours before use.

4.1 Scope of Delivery

Prior to installation, remove the entire delivery scope from the packaging and check that everything is there:

- 1x pressure vessel ①
- 1x pressure vessel lid ③
- 1x filter cartridge ②
- 1x manual


1x carbonate hardness test or total hardness test
1x sticker for Service Pass.
1x reducer 1"-3/4"
1x flush valve with flush hose

If part(s) of the delivery scope are missing, please contact your local SciCan Office (see back of cover).

4.2 Assembly of the pressure vessel and the pressure vessel lid

- Stand on the kick loops ⑨ with both feet.
- Lift the pressure vessel ① and turn it clockwise until the mantle handles ⑪ are over the kick loops ⑨.
- Remove the transport protective cap ⑫ from the filter cartridge.
- Check the O-ring seal of the filter cartridge ② for correct seat in the groove, dirt and damage. **Note:** The cartridge seat has been lubricated with food-safe lubricant at the factory.
- Stand on the kick loops ⑨ with both feet and place the pressure vessel lid ③ on the pressure vessel lid ①. The positioning of the arrow marking on the lid handle ⑬ must line up with the "INSERT" groove.
- Press the pressure vessel lid ③ down and turn clockwise until the lock engages ⑧.

4.3 Assembly of inlet and outlet hoses

Note: The inlet and outlet hoses are not included in the standard scope of delivery. The use of the HYDROS hose sets is recommended.

- Fit inlet hose ⑤ at the inlet of the connector head ④ and outlet hose ⑦ at the outlet of the connector head ④.

Note: Inlet "IN" and outlet "OUT" of the connector head ④ are equipped with O-rings as seals, therefore no additional flat seals may be used here. Make sure that the O-rings are seated correctly.

⚠ Caution: The max. tightening torque at the 1" and 3/4" connections must not exceed 15Nm! Only hose connections with flat seals may be used. Hoses with conical screw connections damage the connections of the filter head and invalidate any guarantee claims! Before assembly, note the direction of flow on the upper side of the filter head, IN = water inlet, OUT = water outlet. Prior to installation, note installation dimensions and operating position (Chapter 10). If no original hoses are used, the 1"- 3/4" ⑭ reducer supplied must be used to ensure correct sealing of the return valve.




5 Commissioning a New Filter

5.1 Identifying the Filter Capacity

- Identify the local carbonate hardness or total hardness using the enclosed carbonate hardness or total hardness test kit.
- Identify the filter capacity using the capacity tables (Chapter 7) taking account of the carbonate or total hardness identified and the bypass setting (Chapter 5.1).

5.2 Flushing and Bleeding the Filter Cartridge

Note: A bucket with a minimum capacity of 10L is needed for flushing/bleeding.

- Place the bucket under the flush hose and open the flush valve (15).
- Open the inlet valve (6) on inlet hose (5) while holding the flush hose in the bucket. Flush with at least 10 litres.
- Close the flush valve (15) and empty the bucket.
- Check system for any leaks.
- Note installation date of the filter system on pressure vessel (1) and next exchange date on the enclosed sticker and attach it to the Service Pass on the pressure vessel (1).

6 Exchanging the Filter Cartridge

The filter cartridge must be replaced after 6–12 months, at the latest 12 months after commissioning, irrespective of the level of exhaustion of the filter system. If the capacity of the filter cartridge has been exhausted (Chapter 7), it must be replaced earlier.

! Caution: During the exchange, carefully examine all dismantled parts! Faulty parts must be exchanged, dirty parts cleaned! Read the operating and safety information (Chapter 3) prior to replacement. After the product has been stored at temperatures below 0 °C, it must be stored with the original packaging open at the ambient temperature of the place of installation for at least 24 hours before use.


6.1 Removing the Filter Cartridge




- Switch off the power supply of the HYDRIM unit (remove plug).
- Close the inlet valve (6) at the inlet hose (5).
- Place the flush hose in a bucket and remove pressure from the filter system by opening the flush valve (15). Collect waste water in the 10L bucket.

Note: If the escaping water is more than 1 litre, the inlet valve (6) is not completely closed or is blocked with scale.

- Stand on the kick loops (9) with both feet while opening the pressure vessel lid (3) by pressing the lock (8) and turning it counter-clockwise as far as it will go.
- Place the pressure vessel lid (3) vertically on the lid handle (13).
- Stand on the kick loops (9) with both feet while turning the pressure vessel (1) counter-clockwise by the mantle handles (11) as far as it will go.
- Take your feet off the kick loops (9) and press the pressure vessel (1) down with both hands on the mantle handles (11).
- Remove used filter cartridge (2) from the pressure vessel (1).
- Place the used filter cartridge (2) in the sink with the connection facing down for drainage (> 5 min).
- Lock the used filter cartridge (2) with the transport protection cap (12) of the new filter cartridge and discard according to the local regulations.

6.2 Inserting the Filter Cartridge

- Check the O-ring seat of the new filter cartridge (2) for correct seat in the groove, dirt and damage. **Note:** The cartridge seat has been lubricated with food-safe lubricant at the factory.
- Check the connector seat of the filter cartridge O-ring (2) in the pressure vessel lid (3) for dirt and damage.
- Place new filter cartridge (2) in the pressure vessel (1).
- Stand on the kick loops (9) with both feet; lift the pressure vessel (1) whole turning clockwise until the mantle handles (11) are over the kick loops (9).
- Stand on the kick loop (9) with both feet and place the pressure vessel lid (3) on the pressure vessel (1). The positioning of the arrow marking on the lid handle (13) must line up with the "INSERT" groove.
- Press the pressure vessel lid (3) down and turn clockwise until the lock engages (8).
- Switch on electrical power supply to the HYDRIM unit.
- Execute the steps described under 5.3 to flush and bleed the new filter cartridge (2).

7 Filter Capacity

Capacity Tables

HYDROS Water Treatment System	
Total hardness in °GH	Volume in litres
4	12.500
5	10.000
6	8.333
7	7.143
8	6.250
9	5.556
10	5.000
11	4.545
12	4.167
13	3.846
14	3.571
15	3.333
16	3.125
17	2.941
18	2.778
19	2.632
20	2.500
21	2.381
23	2.174
25	2.000
28	1.786
31	1.613
35	1.429

Note: The stated capacities are guide values that may vary by $\pm 20\%$ depending on the product volume flow, local water quality and machine type.

8 Repair

Regularly check the filter system for leaks. Regularly check the hoses for kinks. Bent hoses must be replaced.

The complete filter system must be replaced in rotation after a maximum of ten years.
The hoses must be replaced in rotation after a maximum of five years.

! Caution: Prior to exchange read the technical data (Chapter 10) and the operating and safety information (Chapter 3).

Regularly clean the outside of the filter system with a soft, damp cloth.

! Caution: Do not use any substances incompatible with the material (Chapter 3.4) or astringent, abrasive cleaning agents.

9 Troubleshooting

9.1 No water flow

Cause: Water intake closed.

Troubleshooting: Open water intake on the upstream stop valve or inlet valve ⑥ on inlet hose ⑤.

 **Caution:** The following faults may be remedied only by trained and authorised personnel.

9.2 No or low water flow in spite of open water intake

Cause: Mains pressure too low.

Troubleshooting: Check mains pressure.

If the fault continues to occur in spite of adequate mains pressure, check the filter system and filter cartridge and change if necessary.

 **Caution:** Prior to exchange read the technical data (Chapter 10) and the operating and safety information (Chapter 3).

9.3 Leaking screw connections

Cause: Screwed connections not fitted correctly.

Troubleshooting: Check mains pressure. Check all screwed connections and mount according to Chapter 4.

If the fault continues, exchange filter system.

 **Caution:** Prior to exchange read the technical data (Chapter 10) and the operating and safety information (Chapter 3).

9.4 Leak continues after filter has been replaced

Cause: O-ring on filter cartridge not sitting correctly.

Troubleshooting: Check correct seat of the O-ring (Chapter 6.2).

 **Caution:** Prior to dismantling read the data (Chapter 10) and the operating and safety information (Chapter 3).

10 Technical Data

		HYDROS Water Treatment System
Operating pressure		2 bar – max. 6 bar
Operating/water temperature		4°C – 60°C
Ambient temperature during	operation	10°C – 60°C
	storage/transport	–20°C to 60°C
Flow rate with 1 bar pressure loss		850 l/h
Nominal flow according to Norm		300 l/h
Pressure loss at nominal flow		0,45 bar
Weight (dry/wet)		18 kg/24 kg
Dimensions complete system (Width/Depth/Height)		288 mm/255 mm/550 mm
The bending radii of the inlet and outlet hose 2m, DN13, 3/4"–3/4" are 130mm and, depending on the installation orientation and operating space, must be considered in addition to the dimensions of the complete system.		
Operating position		The filter system can be operated either vertically or horizontally.
Inlet connection		G 1"
Outlet connection		G 3/4"

11 Order Numbers

Article	Article Number
HYDROS (complete system with filter cartridge)	01-114436
HYDROS replacement cartridge	01-114450
HYDROS hose set	01-114546
HYDROS Flow meter	01-114547

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Copyright 2025 SciCan Ltd. All rights reserved.