



ENSURE™

Monitoring for Steam Sterilization

Think of it as a routine check up for your sterilizer.

The **Ensure**™ line of Infection Control monitoring products ensures your equipment is performing at its peak. It includes a range of products that are simple, easy to use and help to ensure your patients are being kept safe with every cycle.



Safe

Designed to help keep your patients safe with every cycle.



Easy

Provide the highest level of sterility assurance with easy-to-read indicators, which clearly identify if an instrument load met the expected parameters for sterilization.



Capacity

The 10-well incubator allows you to monitor up to 5 sterilizers at one time, providing you with easy-to-read results that show whether sterility has been achieved or not.



Reliable

Meets the most stringent requirements in today's changing regulatory environment.

Ensure™ Sterilization Monitoring System

Ensure™ Accessories



Ensure™ Chemical Indicators and Integrators

Ensure™ Type 4 and Type 5 indicators are fast and economical solutions for load monitoring. Type 5 Integrators are the most accurate of all chemical indicators and contain a chemical ink which reacts to all three sterilization parameters (time, temperature and pressure).

Ensure™ Sterilization Pouches

Ensure™ Sterilization Pouches have built-in multi-parameter internal/external process indicators that change colour when all three criteria for sterilization have been met - time, temperature and pressure. Once sterilized, contents will remain sterile for up to 3 years provided the integrity of the pouch has not been compromised.



Ensure™ Biological Indicators (BI)

Ensure™ BIs are reliable and self-contained indicators for results in as little as 24 hours. BIs are the only method to ensure effective performance of your sterilizer. The Ensure™ BIs are available in a starter kit system that consists of 25 BIs, the unique incubator, an easy to fill-in log book, and 100 type 5 integrators.



Bowie Dick Test (BD) for Steam Penetration

The Ensure™ BD test is designed for steam vacuum assisted autoclaves to ensure proper steam penetration for a confirmed sterilization process. SciCan's BD test is a unique small form factor test that is easy to use and store. The BD test is available as a starter kit that contains the unique test card holder and 31 test cards.

Part Numbers

Part Number	Description
SCI134	Type 4 Chemical Indicators (250/box)
SCI-T5	Type 5 Chemical Integrators (100/box)
SCI-BI100	Biological Indicators (100/box)
SCI-BISK-120	BI Starter Kit (25 BIs, incubator, log book, 100 type 5 integrators)
SCI-BINC-120	Ensure™ Incubator
SCI-SCP	Sterilization Pouches (200/box)
SCI-BD134	Bowie Dick Test
SCI-BDSK134	Bowie Dick Starter Kit

Manufactured for:
SciCan Ltd.
1440 Don Mills Rd.,
Toronto, ON
Canada, M3B 3P9
Tel: 416-445-1600
Toll free: 1-800-667-7733
customerservice@scican.com

SciCan

www.scican.com