

STATIM G4
Cassette Autoclaves



BORN FAST. BUILT TO LAST.

Back when autoclaves took 45 minutes, STATIM did it in six—making sterilization easy, safe, and predictable. Today, that legacy of speed has grown with an enhanced user interface and digital connectivity, strengthened by 35 years of reliable service. **Solid. Steady. STATIM.**



Trusted by Dentists
in over 100 countries
for 35 years.



Online Monitoring
prevents unscheduled
downtime.



NEW 2-Year Warranty
protects your
investment.



Let's talk Speed

Schedule a consult now!



Capacity

STATIM G4 autoclaves not only offer fast sterilization times but also enough capacity to meet the needs of a busy practice. STATIM 2000 G4 capacity allows for up to 5 pouches and the STATIM 5000 G4 allows for up to 12 pouches.



Gentle

The fresh water feature guarantees fresh water with every cycle. State of the art steam technology is gentle on both solid and hollow instruments.



G4 Technology

G4 equipment can connect and send information to your smart devices so you can have cycle data, instructions, and maintenance notifications at your fingertips.

STATIM 2000 G4 Specifications | G4-121101

Unit size (L x W x H):

19.5" x 16.25" x 6" / 49.5 cm x 41.5 cm x 15 cm

Cassette internal dimensions:

11" x 7" x 1.4" / 28 cm x 18 cm x 3.5 cm

Minimum shelf height for proper leveling:

8.5" / 22 cm

Weight without water:

46 lbs / 21 kg

Power consumption:

115 V, 60 Hz, 11A

	Unwrapped	Wrapped	Rubber / Plastic	Air Drying
Sterilization Time	3.5 minutes	10 minutes	15 minutes	Air drying only
Sterilization Temp	135° C	135° C	121° C	Air drying only
Cycle Time*	6 minutes	14 minutes	20 minutes	Air drying only

*Cycle time is dependent on instrument load and does not include drying.

STATIM 5000 G4 Specifications | G4-201103

Unit size (L x W x H):

23.6" x 16.1" x 7.5" / 60 cm x 41 cm x 19 cm

Cassette internal dimensions:

15" x 7" x 3" / 38 cm x 18 cm x 7.5 cm

Minimum shelf height for proper leveling:

10" / 26 cm

Weight without water:

73 lbs / 33 kg

Power consumption:

115 V, 60 Hz, 11A

	Unwrapped	Wrapped	Rubber / Plastic	Heavy Duty Unwrapped
Sterilization Time	3.5 minutes	6 minutes	35 minutes	6 minutes
Sterilization Temp	132° C	132° C	121° C	132° C
Cycle Time*	9 minutes	17.5 minutes	45 minutes	12 minutes

*Cycle time is dependent on instrument load and does not include drying.



ENSURE

Monitoring for steam sterilization

Easy to use, in-office monitoring for steam sterilization.

SCI-BINC-120



STATMAT

Protect your STATIM

Securely seat cassettes atop your STATIM with a safe and hygienic silicone mat.

01-112867S – For STATIM 2000

01-112868S – For STATIM 5000

Manufactured for:
Dent4You AG
Bahnhofstrasse 2
CH-9435 Heerbrugg

Distributed by:
SciCan Ltd.
1440 Don Mills Road
Toronto, Ontario
M3B 3P9 Canada

Distributed by:
Coltene/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio
44223 US

COLTENEGROUP

SciCan | Your Infection Control Specialist