


Monitoring Instrument Washers



# THINK OF IT AS A ROUTINE CHECK UP

for your washer

The HYDR/M Wash Test Indicators (WTIs) have been designed to provide two-dimensional validation of HYDR/M washing efficacy. It is a simple, easy-to-use indicator that ensures instruments are effectively and consistently cleaned with HIP Cleaning Solution.



#### Safe

Synthetic test soil (non-toxic indicator) designed to monitor the efficacy of HYDR/M washers.



#### Easy

WTI wash off characteristics clearly indicate whether the wash cycle has passed or failed; clean indicator at the end of the cycle indicates an effective cleaning process.



## **Smart Design**

Indicators are folded "back-to-back" and placed inside the holder, monitoring spray performance from two directions.



#### Reliable

Designed to provide verification of a consistently effective cleaning process in HYDR/M washers.

## THE ONLY VALIDATED AND EFFECTIVE **CLEANER IN THE MARKET**

for use in HYDRIM units.

### **HIP Cleaning Solution**

HIP Cleaning Solution is designed to remove any type of protein/debris from instruments. HIP is automatically added at various stages in the cycle to optimize the cleaning and protection of instruments. To maximize HYDR/M cleaning process, it is essential to use the certified and validated HIP Cleaning Solution that works for your HYDR/M unit during the wash cycle.



## **HIP Cleaning Solution for** HYDRIM L110W

Designed for HYDR/M L110W, HIP Cleaning Solution box fits inside L110W, instead of taking up cabinet space. Strong, built-in tube connector prevent leaks.



## HIP Ultra Cleaning Solution for HYDRIM C61W

Designed for HYDR/M C61W, HIP Ultra Cleaning solution pouches fit inside the C61W, instead of taking up cabinet space. Strong, built-in tube connector prevent leaks.

Part Numbers	
Part Number	Description
WTI-SK	Starter Kit (One Holder and 32 Indicators)
WTI-RK	Refill Kit (96 Indicators)
CS-HIPL	HIP Cleaning Solution
CS-HIPC-U	HIP Ultra Solution
C61W-D01	HYDRIM C61W G4
L110W-D01-G4	HYDRIM L110W G4

Manufactured for: SciCan Ltd. 1440 Don Mills Rd., Toronto, ON, Canada, M3B 3P9 Tel: 416-445-1600 Toll free: 1-800-667-7733 customerservice@scican.com

SciCan Inc. 701 Technology Dr., Canonsburg, PA 15317, USA Tel: 724-820-1600 Toll Free: 1-800-572-1211 uscustomerservice@scican.com

Distributed by:

