

AutoFill VistaPure

Look, doc: no hands!

VistaPure's patented AutoFill accessories replenish SciCan STATIM autoclaves with distilled-quality water automatically—enabling truly hands-free convenience.

VistaPure AutoFill accessory kit can be installed easily on any STATIM 2000/G4 or 5000/G4 in just a few minutes!

- ✓ No more purchasing and storing distilled water
- ✓ No more cycle-fault errors caused by forgetting to fill reservoirs
- ✓ No more spills while filling autoclaves manually
- ✓ No more staff time spent filling autoclaves = more time with patients

AutoFill for STATIM (S3020)




VistaPure AutoFill Installation

A 12-foot length of clear polyurethane tubing is already attached to the AutoFill manifold. If the desired distance between the connection to the VistaPure and the autoclave is greater than 12 feet, simply order a longer length of the clear polyurethane tubing, as the distance from the VistaPure to the autoclave is not limited, except by the feed tubing.

Installing AutoFill on a STATIM

The VistaPure AutoFill for STATIM device (S3020) can be installed on any new or older model 2000 or 5000 STATIM autoclave as follows:



Distributed by:
SciCan Ltd.
1440 Don Mills Rd
Toronto, ON M3B 3P9, Canada
customerservice@scican.com

Distributed by:
Coltene/Whaledent Inc.
235 Ascot Pkwy
Cuyahoga Falls, OH 44223, USA

Manufactured by:
Crosstex International, Inc.
1625 Pinnacle Peak Rd
Phoenix, AZ 85027, USA
888-276-7783
HuFriedyGroup.com

www.scican.com

SciCan | Your Infection
Control Specialist