


VistaPure

Water Purification System

Meet your steri-center's new best friend

The VistaPure Water Purification System provides two grades of high-quality water for steam sterilizers, dental bottles, ultrasonic cleaners and the final-rinse cycle of compatible instrument washers. The system connects directly to the city water supply and sends high-quality water anywhere you need it in the steri-center for easy, on-demand use.



Two grades of the water practices need

For steam sterilizers:

Point-and-shoot wand makes filling sterilizers with autoclave-ready water a breeze. No more lifting heavy jugs of water!

For bottle systems:

Mounted faucet provides non-corrosive bottle-quality (hyperfiltered) water for bottle systems.



Instant Status Monitoring

The VistaPure Water Purification System features an inline TDS meter that monitors water quality — and makes it easy to know when it's time to replace your filters.



No Electricity Required


The VistaPure Water Purification System operates on water and air pressure, so no electricity is required to supply all the high-quality water your office needs!



Environmentally Friendly

The VistaPure Water Purification System doesn't require chemicals, and its efficient design minimizes the amount of wastewater produced during filtration.

Optional Accessories



Instrument Washer Adapter Kit
S3010

BRAVO G4 Auxiliary EV Kit
M7C200012


STATCLAVE Adapter Kit
S3015

AutoFill for STATIM Autoclave
S3020

BottleFill Remote Water Dispenser
S3050

V7501
For one autoclave
V7502
For two autoclaves
V7601
For one autoclave

Space Requirements



4.25 Gallon Net Storage Tank
S9111 12" Diameter x 18.25" Height

STANDARD

2 Gallon Net Metal Storage Tank
S9110 8" Diameter x 13" Height

OPTIONAL

9.5 Gallon Net Storage Tank
S9113 15.5" Diameter x 25" Height

OPTIONAL

Manufactured for:
Crosstex International, Inc.
1625 Pinnacle Peak Rd
Phoenix, AZ 85027 USA
T +1 (888) 276-7783
HuFriedyGroup.com

Distributed by:
SciCan Ltd.
1440 Don Mills Rd
Toronto, ON M3B 3P9, Canada
TF +1 (800) 667-7733
T +1 (416) 445-1600
customerservice@scican.com

Distributed by:
Coltene/Whaledent Inc.
235 Ascot Pkwy.
Cuyahoga Falls, OH
44223, USA

For more information, visit www.scican.com

SciCan | Your Infection
Control Specialist

All Vista product names are trademarks of HuFriedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions. HYDRIM and STATIM are registered trademarks, and STATCLAVE and BRAVO are trademarks of SciCan Ltd.